


Lab name: SRI Instruments
 Client: Joann W./Hexion
 Client ID: N10984
 Analysis date: 03/11/2024 11:25:11
 Method: 1ul on-column
 Description: FIDmeth medgain 300C
 Column: 5MXT502 .25mic
 Carrier: H2@2psi=20ml/min
 Integration: Peak sens=95.0 Base sens= 1.0 Min area= 1.00 Standard= 1.000 Sample= 1.000 Tangents=off
 Data file: Lawson497.CHR ()
 Sample: 1000ppm C1-C28 alkanes
 Comments: H2 makeup=0psi air=5psi

Temperature program:

Init temp	Hold	Ramp	Final temp
100.00	0.000	10.000	300.00

Events:

Time	Event
0.000	ZERO
0.000	SOUND



Number	Component	Retention	Area	Internal	Units
14	nC14	1.053	676.7116	0.0000	ppm
16	nC16	2.290	703.8324	0.0000	ppm
18	nC18	3.916	669.2934	0.0000	ppm
22	nC22	5.603	648.2153	0.0000	ppm
24	nC24	7.226	667.1134	0.0000	ppm
26	nC26	8.740	678.7194	0.0000	ppm
28	nC28	10.160	696.5340	0.0000	ppm
7		4740.4195		0.0000	


Lab name: SRI Instruments
 Client: Joann W./Hexion
 Client ID: N10984
 Analysis date: 03/11/2024 11:48:01
 Method: 1ul on-column
 Description: FIDmeth medgain 300C
 Column: 5MXT502 .25mic
 Carrier: H2@2psi=20ml/min
 Integration: Peak sens=95.0 Base sens= 1.0 Min area= 1.00 Standard= 1.000 Sample= 1.000 Tangents=off
 Data file: Lawson498.CHR ()
 Sample: 115mg Caslowax in 40ml CS2
 Comments: H2 makeup=0psi air=5psi

Temperature program:

Init temp	Hold	Ramp	Final temp
100.00	0.000	10.000	300.00

Events:

Time	Event
0.000	ZERO
0.000	SOUND



1ul 115MG CASLOWAX
 DISSOLVED IN 40ML CS2

Number	Component	Retention	Area	Internal	Units
14	nC14	0.940	23.0360	0.0000	ppm
16	nC16	0.000	0.0000	0.0000	ppm
18	nC18	3.696	12.8821	0.0000	ppm
22	nC22	5.423	93.5845	0.0000	ppm
24	nC24	7.043	90.0602	0.0000	ppm
26	nC26	8.550	65.4541	0.0000	ppm
28	nC28	9.960	47.4172	0.0000	ppm
6			332.4341	0.0000	

CASLOWAX
 RETENTION TIME
 COMPARISON
 C16 OVERLAY FOR