



SPECIFICATIONS:

PRODUCT PIN	MR10G-HSSR2-00-0-SRU01
VALVE PIN	MR10G-HSSR2-00-0-R
VALVE TYPE	HIGH PERFORMANCE 5/8GADE. MEDIUM TEMPERATURE
NUMBER OF PORTS	10
PORT SIZE	0.025" (0.64 mm)
PORT CONNECTION	1/16" (1.6 mm)
MAX. WORKING PRESSURE	400 psig (27.6 MPa)
MAX. INITIAL LEAK RATE	8.5E-10 Atm-cm/s
MAX. INITIAL PRESSURE DROP	2 MPa @ 75 acm
STAND-OFF LENGTH	2" (101.6 mm)
ACTUATOR TYPE	ELECTRONIC SINGLE ACTUATION
CONTROL TYPE	DRY CONTACT
ACTUATOR POWER INPUT	24 VDC, 2 A min.
VALVE MAX. WORKING TEMPERATURE	437°F (225°C)
ACTUATOR MAX. WORKING TEMP	140°F (60°C)
MATERIAL	SS-316L VALVE STATOR R2 ROTOR TYPE (MODIFIED FEERX)
TOTAL WEIGHT	2.5 lb (1.1kg)

<p>PROPRIETARY INFORMATION OF AMAPP ON THE INFORMATION CONTAINED HEREIN TRANSMITTED TO ANY THIRD PARTY WITHOUT THE WRITTEN CONSENT OF AMAPP.</p>	
DESCRIPTION	SHI USA 10-Port Relay Valve Kit
FINISH	
SCALE	2:3
SHEET	1 OF 2
REV.	A
DATE	
BY	
CHECKED	
APPROVED	
DATE	

