

MRV10-16R2-0003-02

MRV, 10 Ports,
 Connection: 1/16in Bore: 0.025in,
 Stator: SS Nitronic 60,
 Rotor: R2 Modified PEEK (225°C),
 Customized Mini Electric Actuators: 24 VDC with
 2" Stand-off

Valve S/N: MR11804070
 MEA S/N:
 Date: 24-07-18
 Tested by: VT

Number of Ports: **10**

Connections: 1/16" compression fitting

Bore Size: **.025"** (0.63 mm)

Rotor Material: **R2** – Modified PEEK



Stator Material: SS Nitronic 60

Max. Working Pressure: **80** psig (552 kPa) (5.6 bar)


Maximum Operating Temperature: **225°C** (437°F)

Actuator: Customized MEA


Standoff: 2"



MEA **CW** Position



MEA **CCW** Position



Triggering Connections

Red + Brown: Clockwise

Orange + Brown: Counterclockwise