

SRI Model H2-100 Hydrogen Generator



Make GC quality hydrogen from distilled water.

100ml/minute flow 50 psi outlet pressure

22" x 14" x 20" Shipping weight 52 lbs.

115 or 230VAC power 100 watts

Quiet operation

One Year warranty


The SRI Model H2-100 is a perfect solution for labs which prefer not to have hydrogen cylinders. The H2-100 makes up to 100ml/minute of GC quality hydrogen from distilled water.

The outlet pressure is set to 60 psi but can be adjusted down if required. A quality 1/8" Swagelok style stainless steel outlet fitting makes connection to your GC easy, and a quick shut-off toggle valve provides instant on/off control.

A digital flow display shows H2 usage, and the handy pressure gage confirm proper operation.

To order:

SRI part# 2014-H100 specify voltage (115/230)
100ml/minute hydrogen generator
US\$ 3542.00. **(2022 pricing, prices subject to change, consult most recent price list.)**



SRI Tech Support: 310-214-5092
www.srigc.com

H2100FlyerJan2022