

DETECTOR/TRAPS/CELLS ETC.	SRI US/NON-EUR GC	SCHAMBECK/ORION/SRI EUR
VALVE/MANUAL HEATED VALVES	100W/(1AMP DIODE)	200W/(3AMP DIODE)
P&T TRAPS	100W/(NO DIODE)	200W/(3AMP DIODE)
SPARGE FILTER FOR P&T	100W/(NO DIODE)	100W/(NO DIODE)
SHAKER OR HEADSPACE/VIAL	100W/(1 AMP DIODE)	100W/(1 AMP DIODE)
TCD CELL	100W/(1 AMP DIODE)	200W/(3AMP DIODE)
ECD CELL	100W/(JUMPER)	100W (JUMPER)
RGD CELL	100W/(1 AMP DIODE)	200W/(3 AMP DIODE)
RGD REACTOR	100W/(1 AMP DIODE)	200W/(3 AMP DIODE)
TOGA TRAP	100W/(1 AMP DIODE)	200W/(3 AMP DIODE)
HID SMALL BLOCK	24 Volts / 40 watts	24 Volts / 40 watts
THERMAL DESORBER	150W size 3/8(JUMPER)	150W size 3/8(JUMPER)
5030 P&T auto sampler( soda fountain)	40watts 110/220V	50watts 220V
PTV INJECTOR	200W/(3 AMP DIODE)	200W/(3 AMP DIODE)
HTA INJECTOR	100W/(1 AMP DIODE)	200W/(3 AMP DIODE)
BIO DIESEL DETECTOR	200W/(3 AMP DIODE)	200W/(3 AMP DIODE)
SPLIT SPLITLESS INJECTOR	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
STANDARD POLISHING FILTER	100W (NO DIODE)	100W (NO DIODE)
FID	200W/(3 AMP DIODE)	200W/(3 AMP DIODE)
PID	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
FPD	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
HID	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
TID	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
NPD	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
CCD	N/A	N/A
DELCD	100W (1 AMP DIODE)	200W/(3 AMP DIODE)
PELTIER TRAP	100W (1 AMP DIODE)	100W (1 AMP DIODE)
METHANIZER	200W/(3 AMP DIODE)	200W/(3 AMP DIODE)
510 box	12v 40w	12v 40w