


Model 30 TCD Mini Gas Chromatograph and PeakSimple Data System






Complete GC and Data System
 Super Portable (book size)
 Ambient temperature operation
 No gas cylinders-runs on room air
 Packed or capillary columns
 Includes 2 meter packed column
 TCD detector
 PeakSimple USB Data System
 Made in USA -Two year warranty
Low cost \$3,687.00



The SRI Instruments Model 30 Gas Chromatograph is a simple low-cost system intended for basic GC applications such as measuring Hydrogen in Air or alcohol in water. The extremely small size, USB connection to any Windows computer (laptop) and ability to operate on room air using the included vacuum pump makes the Model 30 perfect for field operations. Helium carrier gas may also be used with the optional carrier gas pressure regulator kit.
(2022 pricing, prices subject to change, consult most recent price list.)

Order:
 0030-1003 Model 30 Basic TCD GC \$3,687.00
 8670-1358 Pressure Regulator w/ Gauge \$513.00



SRI Instruments 310-214-5092
 20720 Earl St. Torrance, CA 90503

email: sales@srigc.com
 website: www.srigc.com